



Features :

- Built with durable steel framework and high-quality ABS & PP plastic
- Three functions controlled by crank (back-rest, leg-rest & hi-low)
- PP head, foot and bed boards are detachable and easy to clean
- 6-column collapsible aluminum alloy protective railing
- 4 pcs universal castor, ø 125mm, individual brake

Technical Specifications:

Dimensions			
Overall dimensions (LxWxH)	2150mmx1050mm x 500-700mm	84.6inx41.3in x 19.7-27.5in	
Hi -Low	700mm-500mm	27.5 -19.7in	
Bed Board dimensions (LxW)	1945 mm x900 mm	76.6inx35.4in	
Diameter of central brake castors	150 mm	5.9 in	
Bed Articulations			
Back -rest		0-60°	
Leg -rest		0-40°	
Trendelenburg and Rev. Trendelenburg		0-10°	
Weight			
Net Weight		139 Kg	306 lbs
Load Capacity		< 250 Kg	<551 lbs
Accessories			
Standard Accessories (Included)		I.V Pole	

Optional Accessories (Not Included)

- Monkey pole
- X-ray radiolucent head and back section