



Features :

- Built with durable steel framework and high-quality ABS & PP plastic
- Three functions easily adjusted by remote control (back-rest, leg-rest and high-low)
- Brake pedal is conveniently located below the foot board

Technical Specifications:

| Dimensions | | |
|-----------------------------------|---------------------------|------------------------------|
| Overall dimensions (LxWxH) | 2160mmx1150mm x 445-655mm | 85.0inx45.3in x 17.5 -25.8in |
| Hi-Low | 655mm-445mm | 25.8-17.5in |
| Bed Board dimensions (LxW) | 1935mmx900mm | 76.2inx35.4in |
| Diameter of central brake castors | 150mm | 5.9in |
| Bed Articulations | | |
| Back -rest | | 0-70° |
| Leg -rest | | 0-35° |
| Weight | | |
| Net Weight | 130 Kg | 287 lbs |
| Load Capacity | <250 Kg | <551 lbs |
| Accessories | | |
| Standard Accessories (Included) | | I.V Pole |

