



Features :

- Built with durable steel framework and high-quality ABS & PP plastic
- Controlled by user-friendly touchpad (replacing traditional keypads, which have limited service life due to mechanical failure)
- Custom bed board improves softness and bending of the bed, increasing the comfort of patients
- Electrically controlled bed positions provide easy adjustment of patient posture, with an indicator light to show the work status
- Back-rest and leg-rest can be adjusted up/down at the same time for more efficient operation
- Quick release lever on back-rest for emergency CPR
- Reset button adjusts the bed to the flat / lowest setting from any position
- Backup battery provides protection for emergency power failure
- Detachable head and foot board made from high-quality PP plastic
- 4 pcs drop-down PP protective railing, controlled by gas spring
- Central brake castors controlled by electric operation
- IV pole mount bracket at all 4 corners of the bed frame



Technical Specifications:

Dimensions			
Overall dimensions (LxWxH)	710mm-530mm	28.0-20.9 in	
Hi -Low	1940mmx900mm	76.4in x 35.4in	
Bed Board dimensions (LxW)	150mm	5.9in	
Diameter of universal central locking castors			
Bed Articulations			
Back -rest		0-65°	
Leg -rest		0-40°	
Trendelenburg and Rev. Trendelenburg		0-10°	
Weight			
Net Weight		148 Kg	326 lbs
Load Capacity		< 250 Kg	<551 lbs
Accessories			
Standard Accessories (Included)		I.V Pole	

Optional Accessories (Not Included)

- Monkey pole
- X-ray radiolucent head and back section

