



Hospital Beds

ToronCare™ Series






EN ISO 13485: 2012 + AC: 2012
FDA QSR 21CFR 820 (cGMP)



Health Canada Establishment Licence



93/42/EEC



Authorized Representatives
issued by
Saudi Food & Drug Authority (SFDA)

About Torontech

The Torontech group is proud to be a Canadian corporation headquartered in Toronto; the commercial, technical, and economic hub of Canada. With offices and an extensive network in the USA as well as other international markets, we are uniquely prepared to provide turn-key solutions for our customers and clients around the globe. At Torontech our quality and variety of products support many industries. Our advantage is sustainable growth in our product lines and providing the best products. We currently offer a wide range of hospital beds, operating tables, birthing beds and surgical lights that are cost-effective, ergonomic, and optimized for all surgical and medical functionality. Also we offer a wide range of Pharmaceutical testers that fully satisfies pharmaceutical functionality as well as cost savings.

Creating comprehensive solutions for our clients has always been a core value of our company. From sales, to order execution, and post-sales support; every staff member is here to help you choose and implement the solution that best suits your unique requirements.

Since inception, we have succeeded in only offering quality manufactured products that are trusted worldwide. We offer the best value for your investment and provide world-class support. In the past few years, Torontech has grown remarkably in numerous industries, but especially in Medical and Pharmaceutical systems. The growth has given us the opportunity to serve many well-known global and North American Organizations, Pharmaceutical, Companies, Hospitals, Clinics & other Health-care establishments.




ToronCare 1076 intelligent bed is our new designed hospital bed to meet multi-function requirements of patients. Controlled by 6pcs of Linak motors with eight functions. The length of bed can extend 20cm.

The embedded operator design is much more convenient for patients to operate. It can be changed to chair position easily.



ToronCare 1078 is a Multifunctional ICU bed with many features like Integrated X-RAY Cassette Holder, In-Bed Scales, Siderail Control, and Supervisor Panel.

This ICU bed is equipped with the automatic compound fallback system which extends the space when positioning the bed in the pelvic and sacral area and significantly reduces pressure on the skin in this zone. At meantime this system will consequently reduces the risk of the formation of pressure ulcers.



ToronCare 1061 is a Multifunctional ICU bed with many features like Integrated X-RAY Cassette Holder, In-Bed Scales, Siderail Control, and Supervisor Panel.

The 1061 Multifunction bed is also equipped with the automatic compound fallback system which extends the space when positioning the bed in the pelvic and sacral area and significantly reduces pressure on the skin in this zone. At meantime this system will consequently reduces the risk of the formation of pressure ulcers.




ToronCare 1060 is an Premium Electric Bed controlled by user-friendly touch pad.

It is built with durable steel framework and high-quality ABS & PP plastic, and electrically controlled bed positions provide easy adjustment of patient posture, with an indicator light to show the work status. The backup battery provides protection for emergency power failure.



ToronCare 1080 - Electric Bed is design for home care with three functions easily adjusted by remote control.

The bed comes with detachable wooden head and foot board, collapsible aluminum protective side rail with user-friendly control switch. It is also equipped with 4 pc stand-board instead of customers, providing increased safety support for patient.



ToronCare 1085 is a detachable Home Care Bed with extra quality and durability.

This bed offers exceptional value for money, but not at the expense of materials or design. It can be assembled, disassembled, and carried by one person without tools in few minutes. Mattresses or anti-decubitus mattress up to a thickness of 8 cm. The 32 cm siderail fall protection height is a great assurance for safety of users.



ToronCare 1010 is a hospital manual bed built with durable steel framework and high-quality ABS & PP plastic. It comes with three functions controlled by crank (back-rest, leg-rest & hi-low). The PP head, foot and bed boards are detachable and easy to clean.



ToronCare 1020 is a hospital manual bed built with durable steel framework and high-quality ABS & PP plastic. It comes with five functions controlled by crank (Back-rest, Leg-rest, Trendelenburg, Rev. Trendelenburg & Hi-Low). The PP head and foot board are detachable and easy to clean.

Hospital Stretchers

Torontech currently offer a wide range if hospital Stretchers that fully satisfies medical functionality as well as cost saving.



Surgical Lights

ToronLight™ is a state of the art surgical light, designed and manufactured according to the highest industry standards. The surgical light provides exceptional site illumination to crate the optimal environment for even the most critical situations.



Operating Tables

At Torontech our quality and variety of products support many industries. Our advantage is sustainable growth in our product lines and providing the best products. We currently offer a wide range of operating tables that fully satisfies surgical and medical functionality as well as cost savings.



Birthing Bed

ToronComfort™ 1000 Birthing Bed has been carefully engineered to provide the most comfortable positions for women during childbirth.





TORONTECH™



Phone: +1 (416) 368-2721



Website: www.Torontech.com



Fax: +1 (416) 981-7652



Email: sales@torontech.com



Address: 7030 Woodbine Ave., Suite #500, Markham, Ontario, L3R 6G2, Canada