

## Features and Applications:

- Manual type transfer stretcher with double hydraulic system.
- Functions are operated by step the foot pedals
- Bright arrow symbols direct the user how to operate clearly.
- Central locking castors are matched with direction castor under the chassis, which helps the user for easy control of stretcher direction.
- Drop-down protective railing, PP bed board produced by blowing molding



## Technical Specifications:

Dimensions	
Overall dimensions (LxW)	1905x750mm (75.0x29.5 in)
Height adjustment (Bed board to ground)	815-580mm (32.0-22.8 in)
Size of bed board (LxW)	1905x610mm (75.0x24.0 in)
Diameter of Central-lock castor	150mm (5.9 in)
Articulations	
Back-rest lifting angle	0-60° ± 5°
Trendelenburg & Rev. Trendelenburg angle	0-13°
Weight	
Net Weight	82 Kg (181 lbs)

## Standard Accessories (Included):

- Mattress
- I.V Pole

