

Features and Applications:

- Manual type transfer stretcher.
- Single crank is used for height adjustment, stretcher can be positioned up and down by turning the crank clockwise and anti-clockwise respectively.
- Bright arrow symbols direct the user how to operate clearly.
- Central locking castors are matched with direction castor under the chassis, which helps the user for easy control of stretcher direction.
- Drop-down protective railing, PP bed board produced by blowing molding



Technical Specifications:

Dimensions	
Overall dimensions (LxW)	2020x740mm (79.5x29.1 in)
Height adjustment (Bed board to ground)	880-580mm (34.6-22.8 in)
Size of bed board (LxW)	1905x610mm (75.0x24.0 in)
Diameter of Central-lock castor	150mm (5.9 in)
Articulation	
Back-rest lifting angle	0-60° ± 5°
Weight	
Net Weight	82 Kg (181 lbs)



Standard Accessories (Included):

- Mattress
- I.V Pole