

Features and Applications:

- Manual connecting Transfer Stretcher, can be used to defend cross-infection.
- Cart A is inside the operation room while Cart B is outside the operating room, both carts are connected by adjusting them to the same height and then pushing Cart B to Cart A, following the orbit part C which can be pushed from one cart to another and vice versa.
- Patient on the mattress can be moved parallel to another bed or flat table by rotating the crank handle on part C.
- Rotary rollover protective railing, insulation bed board.



Technical Specifications:

Dimensions	
Cart A & B connected (LxW)	3880x640mm (152.7x25.2 in)
Height adjustment (Bed board C to ground)	940-690mm (37.0-27.2 in)
Horizontal move range	400mm (15.7 in)
Size of mattress on C (LxW)	1880x555mm (74.0x21.8 in)
Weight	
Net Weight	151 Kg (333 lbs)

Standard Accessories (Included):



Standard Accessories (Included):

- I.V Pole