

INNOVATIVE TECHNOLOGIES

Portable Tablet Hardness Tester

TTPortTAB01



TORONTECH



Portable Tablet Hardness Tester

TTPortTAB01



The TTPortTAB01 Portable Tablet Hardness Tester ensures reliable hardness measurements of tablets, meeting USP requirements. With a load cell capacity of up to 500 N and an optional 1000 N customization, it supports measurements in N, Kgf, and lbf.

This tester features programmable upper (T1) and lower (T2) tolerance limits, displays and prints up to 10 test values, and calculates their average for efficient reporting. It operates on a rechargeable battery and offers RS 232 port connectivity for easy printer integration. Its simple validation process makes it ideal for labs and production environments.

TECHNICAL SPECIFICATIONS

SPECIFICATIONS	DETAILS
Model	TTPortTAB01
Min load	0.1 N
Accuracy	+0.5% full scale
Tablet size	2 mm to 28 mm (length/diameter)
Power	Ni-Hi 8.4V, 1200mA
Dimensions (L x W x H)	355 x 91 x 76 mm
Weight	1.5 Kg (3,3 lbs)

ORDERING INFORMATION

PRODUCTS		
PRODUCT CODE	DESCRIPTION	QUANTITY
0940200	TTPortTAB01 (Body made of S.S.) -230V	1
0930200	TTPortTAB01 (Body made of S.S.) -110V	1


STANDARD ACCESSORIES		
PART NUMBER	DESCRIPTION	QUANTITY
0901A00002	Brush	1
0901A00004	Hook	1

OPTIONAL ACCESSORIES		
PART NUMBER	DESCRIPTION	QUANTITY
0909A00008	40 column printer	1
0901A00001	Brass weight 5 Kgs (Calibrated)	1
090A00011	Starlite thermal paper roll	1




TORONTECH

Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

