

INNOVATIVE TECHNOLOGIES



Benchtop Environmental Chamber

Toron 15L/20L/30L-Series



TORONTECH



Benchtop Environmental Chamber

Toron15L/20L/30L-Series

The Benchtop Environmental Chamber — Toron-Series **provides precise temperature and humidity testing in a compact design**, with capacities of **15L, 30L, and 50L**.

It offers a wide temperature range from **-70°C to +150°C**, with optional humidity control for full climate simulation. Ideal for labs, R&D, and quality control in industries like aviation, automotive, and medical research, it supports stress tests, reliability checks, and weather resistance evaluations on small samples.

Conforms to international standards, including: **IEC60068-2-1, IEC60068-2-2** (Cold & Dry Heat Testing), and **MIL-STD-810F** (High & Low Temperature Testing), this chamber ensures reliable performance under recognized international standards.



FEATURES/ADVANTAGES

► High-Performance Chamber Construction

Built with 1.5mm SUS304 stainless steel interior and cold-rolled steel exterior with electrostatic spray coating for corrosion resistance. TIG-welded seams and drainage holes ensure long-term leak-free use and easy maintenance.

► Precision Air Heating and Circulation

Incorporates stainless steel finned heating tubes, PID-controlled SSR heating, and multi-wing centrifugal fans to deliver fast, stable, and evenly distributed airflow for precise temperature control throughout the chamber.

► Reliable Refrigeration System

Outfitted with a German Secop compressor, sine-wave aluminum fin air-cooled evaporator, and multi-valve electronic expansion throttling. Ensures powerful, efficient cooling and stable low-temperature performance even in demanding environments.

► Accurate Humidity Management

External electric steam humidification using patented generator technology and frozen light dehumidification system. Combined with water quality control, float switch, and micro-vacuum pump for reliable and autonomous humidity regulation.

► Intelligent Control and Monitoring

Equipped with Siemens PLC and a 7-inch TFT touchscreen. Offers real-time data display, program curve settings (up to 20 programs with 99 segments each), alarm diagnostics, data export, and password protection for secure operation.

► Enhanced Observation and Illumination

Dual-coated anti-fog tempered glass viewing window (178–200mm) with integrated transparent heating film. Interior lighting with touchscreen timer provides clear visibility of samples during testing.


► Comprehensive Safety System

Features protection for compressor, fan, heater, and electrical circuits. Includes over-temperature shutdown, alarm recording, visual/audible alerts, and independent secondary temperature protection for maximum user safety.

► Flexible Sample and Cable Access

Includes an adjustable multi-layer sample holder to accommodate various sample sizes and testing heights, along with a Ø50mm cable port sealed with a rubber plug for external device connection without compromising chamber integrity.


TECHNICAL SPECIFICATIONS

SPECIFICATIONS	DETAILS			
Model	Toron15L	Toron30L	Toron50L	
Internal Volume	15L	30L	50L	
Internal Size (mm) (W×D×H)	300 × 200 × 250mm	300 × 300 × 350mm	350 × 350 × 400mm	
External Size (mm) (W×D×H)	550 × 820 × 750mm	550 × 920 × 850mm	600 × 1000 × 1200mm	
Temperature Range	-70°C ~ +150°C			
Temp. Fluctuation	≤±0.3°C			
Temp. Uniformity	≤2.0°C			
Temp. Deviation	≤±2.0°C			
Heating Rate	3°C/min (-40°C ~ 85°C no load, average rate)			
Cooling Rate	2°C/min (85°C ~ -40°C no load, average rate)			
Humidity Range	20% ~ 98%RH 			
Refrigeration Mode	Mechanical compressor refrigeration			
Cooling Method	Air cooling			
Operation Mode	Fixed value control operation, program circulation operation			
Noise	≤72dB (A)			
Remarks	Ambient temperature is +23°C, no-load test, according to international standards.			



TORONTECH

Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652

Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700

Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919

Email: sales@torontech.com



USA | Canada | EU | UAE | GCC | India | Africa | Latin America | Asia

TORONTECH.com