



Features:

- High quality poly surface for decubitus prevention, removable, washable, impermeable
- Kidney bridge up
- Backlit display indicating status information, e.g. currently chosen position, and status of battery charging, it also ensures a safe and easy handling even in the dark operation room.
- Additional control at the bottom also is with backlit
- Stainless steel inlay plates
- Diam 125mm castors with locking, universal castors
- Many selections of accessories, extensions, and positioning aids

Technical Data

Total length (with additional footplate)	2200mm
Total width:	600 mm
Net Width:	555mm
Length of Head section:	300mm
Length of Back section:	380mm
Length of Main section:	1080mm
Length of Leg section:	580mm
Supporting bottom:	910x625x140mm
Height adjustment range:	700-1000mm
Kidney support up:	120mm
Movement backward to forward:	350mm
Back section up:	70 degree
Back section down:	22 degree
Trendelenburg:	30 degree
Anti-trendelenburg:	-30 degree
Lateral adjustment:	22 degree
Leg section inclination:	+15degree till -90degree
Spread angle of leg plates:	100 degree
Head section inclination:	+45degree till -90degree
Net weight of overall table:	280 kgs
Patient horizontal max load:	300 kgs
Max load weight:	320 kgs