



## Advantages

### Quality and Durability

The homecare bed is sturdily constructed to match with the rigours of daily care. High-quality components make this a particularly durable, sturdy and low maintenance care bed that reduces further costs.

### Cost-Effective

This bed offers exceptional value for money, but not at the expense of materials or design.

### Easy for Assembling & Carrying

The homecare bed can be assembled, disassembled, and carried by one person without tools in few minutes.

### Safety Side Fall Protection

Mattresses or anti-decubitus mattress up to a thickness of 8 cm. The 32 cm siderail fall protection height is a great assurance for safety of users.

### Safety Side Fall Protection

Functions	
Backrest Angle:	0-72°
Knee Section Angle:	0-40°
Height Adjust:	335 -725mm
Trendelenburg and Reverse Trendelenburg angle	0~15°
Configurations	
Size:	2180*1000* (335 -725 ) mm
Bed Board:	Wooden bar with high loading weight
Head and Foot Board:	Wooden
Guardrail:	Aluminum with surface processing, which is durable and warping preventing.
Castor:	Universal castors with brakes
Optional Configuration	
Grip pole, Movable table, Bedside cabinet, Mattress	
Warranty	
Motor: 2 Years, Other Parts:1 year	
Specialty	
The bed is foldable which saves space for shipping and storage.	
It is easy to install the beds by one person without the help of any tools, which saves time and energy.	

