



### Features :

- Built with durable steel framework and high-quality ABS & PP plastic
- Five functions controlled by crank (back-rest, leg-rest, Trendelenburg, Rev. Trendelenburg & hi-low)
- PP head and foot board are detachable and easy to clean
- 4 pcs drop-down PP protective railing, controlled by gas spring
- Brake pedal is conveniently located below the foot board

### Technical Specifications:

Models	ToronCare 1020™	
	Dimensions	
Overall dimensions (LxWxH)	2200mmx1050mm x 590-680mm 86.6inx41.3in x 23.2-26.7in	
Hi-Low	680mm-590mm 26.7-23.2in	
Bed Board dimensions (LxW)	1950mmx900mm 76.7inx35.4in	
Diameter of central brake castors	150mm 5.9in	
	Bed Articulations	
Back -rest	0-70°	
Leg -rest	0-40°	
Trendelenburg and Rev. Trendelenburg	0-10°	
	Weight	
Net Weight	117 Kg 258 lbs	
	Accessories	
Standard Accessories (Included)	I.V Pole	

